# NOTICE

## Masque Thermo-Thérapeutique *Hot Cold Therapeutic Gel Mask*

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*The thermo-therapeutic gel mask of Sun & Sia offers you a hot or cold care at your fingertips.*

*Our mask can be used to calm and comfort irritated skin from sunburn or aesthetic invasive treatments. Relieve stress related tension and migraines, sinus pressure. It can help to reduce puffy eyes and dark circles, also helps to lift the face and prevent the facial skin from loosening.*

***How to use:***

***For cold use:*** place the facial mask in the refrigerator for 2 hours before use. Fit the facial mask around the head, placing the gel portion over the eyes and face, adjust as needed for best fit.

***For hot use:*** Place in hot water, 80°C max for 4-6 minutes. Check the temperature of the mask before use. Fit the facial mask around the head, placing the gel portion over the eyes and face, adjust as needed for best fit.

***Benefits:***

- ✓ *Calm & comfort irritated skin (sunburn, aesthetic invasive treatments).*
- ✓ *Relieve stress related tension and migraines, sinus pressure.*
- ✓ *Helps to reduce puffy eyes and dark circles.*
- ✓ *Lift the face and prevent skin from loosening.*
- ✓ *Refine pores and make the complexion more radiant.*

***Caution:***

*For external use only.*

*Discard mask if damaged or leaking.*

*If contents of the mask makes contact with eyes or skin, wash thoroughly with water.*

*Do not use on infants.*

*Keep out of reach of children.*

---